

Anti-Slc34A2 Polyclonal Antibody For Wb And Elisa - O95436

Item Number: CM0026088

Introduction

SLC34A2 Rabbit Polyclonal Antibody validated for Western blot (WB) and ELISA, reacting with Human, Mouse, and Rat. Suitable for investigating sodium-dependent phosphate transport in lung alveolar type II cells.

[Learn More](#)

Parameter	Value
Product Name	SLC34A2 Rabbit pAb
Remarks/Alias	PULAM; NPT1Ib; NaPi2b; NAPI-3B; NAPI-I1b; SLC34A2
Species	Human
GeneID (Human)	10568
GeneID	10568
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 234-362 of human SLC34A2 (NP_001171469.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95436
Protein Weight	76kDa
Shipping	Ice bag