

Anti-Slc34A3 Rabbit Polyclonal Antibody For Wb/Elisa - Q8N130

Item Number: CM0018151

Introduction

Rabbit polyclonal antibody against human SLC34A3 (NaPi-2c, NPT2C) validated for Western blot and ELISA. Reacts with human, mouse, and rat kidney samples. Useful for renal phosphate transport and HHRH research.

[Learn More](#)

Parameter	Specification
Product Name	SLC34A3 Rabbit pAb
Remarks/Alias	HHRH; NPT2C; NPT1c; SLC34A3
Species	Human
GeneID (Human)	142680
GeneID	142680
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-117 of human SLC34A3 (NP_543153.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N130
Protein Weight	64kDa
Shipping	Ice bag