

Anti-Slc35A1 Rabbit Polyclonal Antibody For Wb, Elisa - P78382

Item Number: CM0014369

Introduction

Rabbit polyclonal antibody targeting human SLC35A1 (CMP-sialic acid transporter). Validated for WB and ELISA; cross-reacts with human and mouse. Suitable for glycosylation and Golgi transport research.

[Learn More](#)

Specification	Details
Product Name	SLC35A1 Rabbit pAb
Remarks/Alias	CST; hCST; CDG2F; CMPST; SLC35A1
Species	Human
Gene ID (Human)	10559
Gene ID	10559
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human SLC35A1 (NP_006407.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P78382
Protein Weight	37kDa
Shipping	Ice bag