

Anti-Slc35A2 Polyclonal Antibody For If/Icc And Elisa - P78381

Item Number: CM0029066

Introduction

Rabbit polyclonal antibody against human SLC35A2 (UDP-galactose translocator).
Validated for IF/ICC and ELISA. Suitable for studying glycosylation and Golgi transport. - P78381.

[Learn More](#)

Parameter	Value
Product Name	SLC35A2 Rabbit pAb
Remarks/Alias	UGT; CDGX; UGAT; UGT1; UGT2; UGTL; CDG2M; UGALT; UDP-Gal-Tr; SLC35A2
Species	Human
GeneID	7355
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-242 of human SLC35A2 (NP_001027460.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P78381
Protein Weight	41kDa
Shipping	Ice bag