

Anti-Slc35D3 Polyclonal Antibody For Wb, Elisa - Q5M8T2

Item Number: CM0023235

Introduction

SLC35D3 rabbit polyclonal antibody (Gene ID: 340146) validated for WB and ELISA. Detects human SLC35D3, a probable UDP-glucose transporter involved in synaptic vesicle function and autophagy.

[Learn More](#)

Parameter	Value
Product Name	SLC35D3 Rabbit pAb
Remarks/Alias	FRCL1; SLC35D3
Species	Human
Gene ID (Human)	340146
Gene ID	340146
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 310-390 of human SLC35D3 (NP_001008783.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q5M8T2
Protein Weight	44kDa
Shipping	Ice bag