


Anti-Slc37A4 Polyclonal Antibody For Wb, Ihc-P, Elisa - 043826

Item Number: CM0017048



Introduction

Rabbit polyclonal antibody targeting human SLC37A4 (G6PT) for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Supports glucose homeostasis and glycogen storage disease research.

[Learn More](#)

Parameter	Value
Product Name	SLC37A4 Rabbit pAb
Remarks/Alias	G6PT; SPX4; CDG2W; G6PT1; G6PT2; G6PT3; GSD1b; GSD1c; GSD1d; TRG19; TRG-19; PRO0685; SLC37A4
Species	Human
GeneID (Human)	2542
GeneID	2542
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 330-429 of human SLC37A4 (NP_001458.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43826
Protein Weight	46kDa
Shipping	Ice bag