

Anti-Slc39A11 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Bwy7

Item Number: CM0021302

Introduction

High-quality rabbit polyclonal antibody against mouse Slc39a11 (ZIP-11), a zinc transporter. Validated for western blot, immunohistochemistry, and ELISA. Cross-reacts with human and rat. Ideal for metal ion homeostasis studies.

[Learn More](#)

Parameter	Value
Product Name	Slc39a11 Rabbit pAb
Remarks/Alias	ZIP-11; 1810074D23Rik; Slc39a11
Species	Mouse
GeneID	69806
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 101-174 of mouse Slc39a11 (NP_001159975.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8BWY7
Protein Weight	35kDa
Shipping	Ice bag