

Anti-Slc39A2 Polyclonal Antibody For Wb, Elisa - Q9Np94

Item Number: CM0021225

Introduction

Rabbit polyclonal antibody targeting human SLC39A2 (ZIP2), validated for WB and ELISA with cross-reactivity to rat. Suitable for zinc transport studies and keratinocyte differentiation research.

[Learn More](#)

Parameter	Value
Product Name	SLC39A2 Rabbit pAb
Remarks/Alias	6A1; ZIP2; ETI-1; ZIP-2; SLC39A2
Species	Human
GeneID (Human)	29986
GeneID	29986
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SLC39A2 (NP_055394.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q9NP94
Protein Weight	33kDa
Shipping	Ice bag