

Anti-Slc39A4 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q6P5W5

Item Number: CM0028349

Introduction

Rabbit polyclonal antibody targeting human SLC39A4 (ZIP4), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Zinc transporter essential for dietary zinc uptake; ideal for zinc metabolism research.

[Learn More](#)

Parameter	Value
Product Name	SLC39A4 Rabbit pAb
Remarks/Alias	AEZ; ZIP4; AWMS2; SLC39A4
Species	Human
GeneID	55630
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-327 of human SLC39A4 (NP_570901.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6P5W5
Protein Weight	68kDa
Shipping	Ice bag