

Anti-Slc39A6 Polyclonal Antibody For Wb, Elisa - Q13433

Item Number: CM0025608

Introduction

Rabbit polyclonal antibody targeting human SLC39A6 (ZIP6/LIV-1), a zinc transporter involved in EMT. Validated for Western blot (85 kDa) and ELISA. Cross-reacts with mouse and rat. Ideal for cancer, immunology, and zinc homeostasis research.

[Learn More](#)

Parameter	Value
Product Name	SLC39A6 Rabbit pAb
Remarks/Alias	LIV1; ZIP6; LIV-1; SLC39A6
Species	Human
GeneID	25800
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 170-320 of human SLC39A6 (NP_001092876.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13433
Protein Weight	85kDa
Shipping	Ice bag