

Anti-Slc39A7 (Zip7) Rabbit Polyclonal Antibody For Wb And Elisa - Q92504

Item Number: CM0029798



Introduction

Rabbit polyclonal antibody targeting SLC39A7 (ZIP7), an ER/Golgi zinc transporter critical for zinc homeostasis, B cell development, and cell signaling. Validated for WB and ELISA in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	SLC39A7 Rabbit pAb
Remarks/Alias	KE4; AGM9; HKE4; ZIP7; RING5; H2-KE4; D6S115E; D6S2244E; SLC39A7
Species	Human
Gene ID (Human)	7922
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 235-380 of human SLC39A7 (NP_008910.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92504
Protein Weight	50kDa
Shipping	Ice bag