

Anti-Slc39A8 Rabbit Polyclonal Antibody For Wb, Elisa - Q9C0K1

Item Number: CM0024087

Introduction

Rabbit polyclonal antibody targeting human SLC39A8 (ZIP8), validated for Western blot and ELISA with cross-reactivity to rat. Suitable for studying zinc/manganese transport, immune regulation, and metal homeostasis.

[Learn More](#)

Parameter	Value
Product Name	SLC39A8 Rabbit pAb
Remarks/Alias	ZIP8; CDG2N; PP3105; BIGM103; LZT-Hs6; SLC39A8
Species	Human
GeneID (Human)	64116
GeneID	64116
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human SLC39A8 (NP_071437.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9C0K1
Protein Weight	50kDa
Shipping	Ice bag