


Anti-Slc3A1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q07837

Item Number: CM0028249



Introduction

High-quality rabbit polyclonal antibody against human SLC3A1, validated for WB, IHC-P, and ELISA. Detects endogenous levels in human, mouse, and rat. Ideal for amino acid transport and kidney disease research.

[Learn More](#)

Parameter	Value
Product Name	SLC3A1 Rabbit pAb
Remarks/Alias	D2H; ATR1; NBAT; RBAT; CSNU1; SLC3A1
Species	Human
Gene ID (Human)	6519
Gene ID	6519
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 416-685 of human SLC3A1 (NP_000332.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07837
Protein Weight	79kDa
Shipping	Ice bag