

Anti-Slc3A2/Cd98Hc Rabbit Polyclonal Antibody For Wb, Elisa - P08195

Item Number: CM0024124



Introduction

Rabbit polyclonal antibody targeting SLC3A2/CD98hc, suitable for Western blot (WB) and ELISA. Cross-reacts with human, mouse, and rat. Ideal for amino acid transport, integrin signaling, and viral entry research.

[Learn More](#)

Product Name	SLC3A2/CD98hc Rabbit pAb
Remarks/Alias	4F2; CD98; MDU1; 4F2HC; 4T2HC; NACAE; CD98HC
Species	Human
Gene ID (Human)	6520
Gene ID	6520
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 461-630 of human SLC3A2/CD98hc (NP_002385.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08195
Protein Weight	58kDa
Shipping	Ice bag