

Anti-Slc3A2/Cd98Hc Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P08195

Item Number: CM0027312



Introduction

Rabbit polyclonal antibody against SLC3A2/CD98hc, validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Functions as amino acid transporter chaperone and is involved in integrin signaling and viral entry.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] SLC3A2/CD98hc Rabbit pAb
Remarks/Alias	4F2; CD98; MDU1; 4F2HC; 4T2HC; NACAE; CD98HC; SLC3A2/CD98hc
Species	Human
GeneID (Human)	6520
GeneID	6520
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 206-630 of human SLC3A2/CD98hc (NP_002385.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08195
Protein Weight	58kDa
Shipping	Ice bag