


Anti-Slc40A1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Np59

Item Number: CM0013742



Introduction

Specific anti-SLC40A1 rabbit polyclonal antibody for Western blot and ELISA. Detects human and mouse ferroportin (Q9NP59), a key iron exporter. Suitable for IVD research and assay development.

[Learn More](#)

Parameter	Value
Product Name	SLC40A1 Rabbit pAb
Remarks/Alias	FPN; FPN1; HFE4; MTP1; IREG1; MST079; MSTP079; SLC11A3; SLC40A1
Species	Human
GeneID (human)	30061
GeneID	30061
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 250-300 of human SLC40A1 (NP_055400.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NP59
Protein Weight	63kDa
Shipping	Ice bag