

Anti-Slc44A1 Polyclonal Antibody For Wb, Elisa - Q8Wwi5

Item Number: CM0017695



Introduction

Rabbit polyclonal antibody targeting human SLC44A1 (choline transporter-like protein 1). Validated for WB and ELISA, cross-reacts with mouse and rat. Suitable for choline transport and membrane synthesis studies.

[Learn More](#)

Parameter	Value
Product Name	SLC44A1 Rabbit pAb
Remarks/Alias	CD92; CTL1; CDW92; CHTL1; CONATOC; SLC44A1
Species	Human
GeneID (human)	23446
GeneID	23446
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 55-210 of human SLC44A1 (NP_536856.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WWI5
Protein Weight	73kDa
Shipping	Ice bag