

Anti-Slc46A1 Rabbit Polyclonal Antibody For Wb And Elisa - Q96Nt5

Item Number: CM0030125

Introduction

Rabbit polyclonal antibody raised against N-terminal region of human SLC46A1 (PCFT/HCP1). Suitable for WB and ELISA applications. Detects endogenous SLC46A1 in human, mouse, and rat. For folate transport and heme trafficking research.

[Learn More](#)

Parameter	Value
Product Name	SLC46A1 Rabbit pAb
Remarks/Alias	G21; HCP1; PCFT; hPCFT; HsPCFT; SLC46A1
Species	Human
GeneID (Human)	113235
GeneID	113235
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human SLC46A1 (NP_001229295.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96NT5
Protein Weight	50kDa
Shipping	Ice bag