


Anti-Slc47A1 (Mate1) Rabbit Polyclonal Antibody For Wb And Elisa - Q96FI8

Item Number: CM0017336



Introduction

Human SLC47A1/MATE1 rabbit polyclonal antibody validated for WB and ELISA. Recognizes human and rat MATE1 multidrug efflux transporter; ideal for drug transport, excretion, and toxicology research.

[Learn More](#)

Parameter	Value
Product Name	SLC47A1 Rabbit pAb
Remarks/Alias	MATE1; SLC47A1
Species	Human
Gene ID (Human)	55244
Gene ID	55244
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 460-550 of human SLC47A1 (NP_060712.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q96FL8
Protein Weight	62kDa
Shipping	Ice bag