

Anti-Slc4A1 Polyclonal Antibody For Wb, If-P, Elisa - P02730

Item Number: CM0019078



Introduction

Rabbit polyclonal antibody targeting the Band 3 anion transport protein (SLC4A1/AE1). Validated for WB, IF-P, and ELISA applications; cross-reacts with human, mouse, and rat. Ideal for erythrocyte membrane, renal physiology, and malaria invasion studies.

[Learn More](#)

Parameter	Value
Product Name	SLC4A1 Rabbit pAb
Remarks/Alias	DI; FR; SW; WD; WR; AE1; CHC; SAO; WD1; BND3; EPB3; SPH4; CD233; EMPB3; RTA1A; SLC4A1
Species	Human
Gene ID (Human)	6521
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-353 of human SLC4A1 (NP_000333.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02730
Protein Weight	102kDa
Shipping	Ice bag