


Anti-Slc4A11 Polyclonal Antibody For Wb, Elisa - Q8Nbs3

Item Number: CM0019644



Introduction

Rabbit polyclonal antibody against human SLC4A11 (NaBC1) for WB, ELISA. Cross-reacts with human and mouse. Useful for research on boron transport, corneal dystrophy, and ion homeostasis.

[Learn More](#)

Parameter	Value
Product Name	SLC4A11 Rabbit pAb
Remarks/Alias	BTR1; CHED; CDPD1; CHED2; NABC1; dj79416.2; SLC4A11
Species	Human
GeneID (Human)	83959
GeneID	83959
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human SLC4A11 (NP_001167560.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8NBS3
Protein Weight	98kDa
Shipping	Ice bag