


Anti-Slc4A2 Polyclonal Antibody For Wb, Elisa - P04920

Item Number: CM0029277



Introduction

Rabbit polyclonal antibody targeting human anion exchange protein 2 (SLC4A2/AE2). Validated for WB and ELISA, this reagent detects a 137kDa chloride/bicarbonate antiporter critical for osteoclast function and intracellular pH regulation.

[Learn More](#)

Parameter	Value
Product Name	SLC4A2 Rabbit pAb
Remarks/Alias	AE2; HKB3; BND3L; NBND3; EPB3L1; SLC4A2
Species	Human
Gene ID (Human)	6522
Gene ID	6522
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P04920
Protein Weight	137kDa
Shipping	Ice bag