

Anti-Slc4A4 / Nbc Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y6R1

Item Number: CM0029673



Introduction

Rabbit polyclonal antibody against human SLC4A4/NBC (electrogenic sodium bicarbonate cotransporter 1), validated for WB, IF/ICC, ELISA; reacts with human, mouse, rat. Ideal for ion transport and pH regulation studies.

[Learn More](#)

Parameter	Specification
Product Name	SLC4A4 / NBC Rabbit pAb
Remarks/Alias	KNBC; NBC1; NBC2; pNBC; HNBC1; NBCe1; hhNMC; kNBC1; SLC4A5; NBCe1-A; SLC4A4 / NBC
Species	Human
GeneID (Human)	8671
GeneID	8671
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 978-1078 of human SLC4A4 / NBC (NP_001091954.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6R1
Protein Weight	121kDa
Shipping	Ice bag