

Anti-Slc4A5 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9By07

Item Number: CM0016357

Introduction

SLC4A5 Rabbit polyclonal antibody validated for WB, IHC-P, ELISA. Detects human, mouse, rat electrogenic sodium bicarbonate cotransporter 4 (NBCe2) for transport research.

[Learn More](#)

Parameter	Value
Product Name	SLC4A5 Rabbit pAb
Remarks/Alias	NBC4; NBCe2; SLC4A5
Species	Human
Gene ID (Human)	57835
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1042-1121 of human SLC4A5 (NP_597812.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BY07
Protein Weight	126kDa
Shipping	Ice bag