


Anti-Slc4A7 Polyclonal Antibody For Wb - Q9Y6M7

Item Number: CM0016002



Introduction

CamelBio's SLC4A7 Rabbit Polyclonal Antibody (WB, ELISA) detects human and mouse SLC4A7. Ideal for studying sodium bicarbonate cotransport and pH regulation in research.

[Learn More](#)

Parameter	Value
Product Name	SLC4A7 Rabbit pAb
Remarks/Alias	NBC2; NBC3; SBC2; NBCN1; SLC4A6; SLC4A7
Species	Human
Gene ID (Human)	9497
Gene ID	9497
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1151-1214 of human SLC4A7 (NP_003606.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y6M7
Protein Weight	136kDa
Shipping	Ice bag