

Anti-Slc4A8 Rabbit Polyclonal Antibody For Wb, Elisa - Q2Y0W8

Item Number: CM0018716



Introduction

Rabbit polyclonal antibody targeting human SLC4A8, suitable for WB and ELISA. Cross-reacts with human, mouse, and rat. Useful for pH regulation and ion transport research.

[Learn More](#)

Parameter	Value
Product Name	SLC4A8 Rabbit pAb
Remarks/Alias	NBC3, NDCBE, SLC4A8
Species	Human
Gene ID (Human)	9498
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-90 of human SLC4A8 (NP_001035049.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q2Y0W8
Protein Weight	123kDa
Shipping	Ice bag