

Anti-Slc5A1 Rabbit Polyclonal Antibody For Wb, Elisa - P13866

Item Number: CM0014893



Introduction

SLC5A1 rabbit polyclonal antibody validated for WB and ELISA, specifically reactive with human samples. Ideal for investigating sodium/glucose cotransporter 1 in intestinal glucose absorption and reproductive biology.

[Learn More](#)

Parameter	Value
Product Name	SLC5A1 Rabbit pAb
Remarks/Alias	NAGT; SGLT1; SGLT-1; D22S675; SLC5A1
Species	Human
GeneID	6523
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 400-500 of human SLC5A1 (NP_000334.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P13866
Protein Weight	73kDa
Shipping	Ice bag