

Anti-Slc5A2 Rabbit Polyclonal Antibody For Wb, If, Ihc, Elisa - P31639

Item Number: CM0022832



Introduction

Rabbit polyclonal antibody against human SLC5A2 (SGLT2), validated for WB, IF, IHC, and ELISA. Detects endogenous levels in human, mouse, and rat samples, with a calculated molecular weight of 73 kDa.

[Learn More](#)

Parameter	Description
Product Name	SLC5A2 Rabbit pAb
Remarks/Alias	SGLT2; SLC5A2
Species	Human
GeneID (Human)	6524
GeneID	6524
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 564-624 of human SLC5A2 (NP_003032.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31639
Protein Weight	73kDa
Shipping	Ice bag