

Anti-Slc5A3 Polyclonal Antibody For Wb, Elisa - P53794

Item Number: CM0019021

Introduction

Rabbit polyclonal antibody against SLC5A3, validated for WB and ELISA. Reacts with human, mouse, rat. Suitable for studying sodium/myo-inositol cotransport, osmotic stress, and ion channel regulation. Predicted weight 80kDa.

[Learn More](#)

Parameter	Value
Product Name	SLC5A3 Rabbit pAb
Remarks/Alias	BCW2; SMIT; SMIT1; SMIT2; SLC5A3
Species	Human
Gene ID (Human)	6526
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 535-695 of human SLC5A3 (NP_008864.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P53794
Protein Weight	80kDa
Shipping	Ice bag