

Anti-Slc5A5 Polyclonal Antibody For Wb, Elisa - Q92911

Item Number: CM0027400

Introduction

Rabbit polyclonal antibody raised against human SLC5A5 (sodium/iodide cotransporter). Suitable for Western blot and ELISA, detecting human, mouse, and rat NIS. Ideal for thyroid research.

[Learn More](#)

Parameter	Value
Product Name	SLC5A5 Rabbit pAb
Remarks/Alias	NIS; TDH1; SLC5A5
Species	Human
Gene ID (Human)	6528
Gene ID	6528
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 544-643 of human SLC5A5 (NP_000444.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92911
Protein Weight	69kDa
Shipping	Ice bag