


Anti-Slc5A7 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q9Gzv3

Item Number: CM0028400



Introduction

High-affinity choline transporter SLC5A7 (CHT1) rabbit polyclonal antibody validated for WB, IF-P, ELISA. Suitable for human, mouse, rat samples. Ideal for cholinergic neuron research.

[Learn More](#)

Parameter	Value
Product Name	SLC5A7 Rabbit pAb
Remarks/Alias	CHT; CHT1; CMS20; HMN7A; DHMNVP; SLC5A7
Species	Human
GeneID (Human)	60482
GeneID	60482
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 503-580 of human SLC5A7 (NP_068587.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZV3
Protein Weight	63kDa
Shipping	Ice bag