

Anti-Slc5A9 Rabbit Polyclonal Antibody For Wb And Elisa - Q2M3M2

Item Number: CM0023368

Introduction

Anti-SLC5A9 rabbit polyclonal antibody validated for WB and ELISA. Reacts with human and mouse SLC5A9/SGLT4 (74 kDa), a sodium/glucose cotransporter involved in mannose transport. Suitable for ion transport and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	SLC5A9 Rabbit pAb
Remarks/Alias	SGLT4; SLC5A9
Species	Human
GeneID (Human)	200010
GeneID	200010
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 230-320 of human SLC5A9 (NP_001011547.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q2M3M2
Protein Weight	74kDa
Shipping	Ice bag