


Anti-Slc6A1 Polyclonal Antibody For Wb, Elisa - P30531

Item Number: CM0018037



Introduction

Rabbit polyclonal antibody against human SLC6A1 (GAT1), raised against a synthetic peptide within amino acids 500-599. Validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Ideal for neuroscience and GABA transporter research.

[Learn More](#)

Parameter	Value
Product Name	SLC6A1 Rabbit pAb
Remarks/Alias	MAE; GAT1; GABATR; hGAT-1; GABATHG; SLC6A1
Species	Human
Gene ID (Human)	6529
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 500-599 of human SLC6A1 (NP_003033.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30531
Protein Weight	67kDa
Shipping	Ice bag