

Anti-Slc6A14 Polyclonal Antibody For Wb, Elisa - Q9Un76

Item Number: CM0015329

Introduction

Rabbit polyclonal antibody against human SLC6A14, validated for WB and ELISA. Cross-reacts with mouse. The target is a sodium- and chloride-dependent neutral and basic amino acid transporter B(0+).

[Learn More](#)

Parameter	Value
Product Name	SLC6A14 Rabbit pAb
Remarks/Alias	BMIQ11; SLC6A14
Species	Human
Gene ID (Human)	11254
Gene ID	11254
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 131-234 of human SLC6A14 (NP_009162.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UN76
Protein Weight	72kDa
Shipping	Ice bag