

Anti-Slc6A2 (Norepinephrine Transporter) Rabbit Polyclonal Antibody For Wb And Elisa - P23975

Item Number: CM0012951



Introduction

Rabbit polyclonal antibody against human SLC6A2 (Norepinephrine transporter, NET) validated for Western blotting and ELISA. Cross-reacts with human, mouse, and rat. Recombinant immunogen. Suitable for neuroscience research and transporter studies.

[Learn More](#)

Parameter	Value
Product Name	SLC6A2 Rabbit pAb
Remarks/Alias	NET; NAT1; NET1; SLC6A5; SLC6A2
Species	Human
Gene ID (Human)	6530
Gene ID	6530
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-235 of human SLC6A2 (NP_478066.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23975
Protein Weight	69kDa
Shipping	Ice bag