

Anti-Slc6A2 Rabbit Monoclonal Antibody For Wb, Elisa - P23975

Item Number: CM0006615



Introduction

High-specificity rabbit monoclonal antibody targeting human SLC6A2 (NET, norepinephrine transporter) validated for WB and ELISA applications. Ideal for neurotransmitter reuptake studies and sympathetic nervous system research.

[Learn More](#)

Specification	Details
Product Name	SLC6A2 Rabbit mAb
Remarks/Alias	NET; NAT1; NET1; SLC6A5; SLC6A2
Species	Human
GeneID (human)	6530
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human SLC6A2 (NP_001034.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P23975
Protein Weight	69kDa
Shipping	Ice bag