

Anti-Slc6A5 (Glyt2) Polyclonal Antibody For Wb And Elisa - Q9Y345

Item Number: CM0012891

Introduction

Rabbit polyclonal antibody targeting human SLC6A5 glycine transporter 2. Validated for Western blotting and ELISA, detects endogenous GlyT2 in human, mouse, and rat samples. Suitable for neuroscience and transporter research.

[Learn More](#)

Parameter	Value
Product Name	SLC6A5 Rabbit pAb
Remarks/Alias	NET1, GLYT2, HKPX3, GLYT-2, SLC6A5
Target Species	Human
GeneID	9152 (Human)
Immunogen	Recombinant fusion protein containing amino acids 1-200 of human SLC6A5 (NP_004202.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y345
Protein Weight	87 kDa
Shipping	Ice bag