


Anti-Slc6A8 Rabbit Polyclonal Antibody For Wb And Elisa - P48029

Item Number: CM0022954



Introduction

Rabbit polyclonal antibody targeting human SLC6A8 (creatine transporter), suitable for WB and ELISA; cross-reacts with mouse. Supports research on creatine transport in brain and muscle.

[Learn More](#)

Parameter	Value
Product Name	SLC6A8 Rabbit pAb
Remarks/Alias	CRT; CT1; CRT1; CRTR; CTR5; CCDS1; CRT-1; SLC6A8
Species	Human
GeneID (Human)	6535
GeneID	6535
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 539-635 of human SLC6A8 (NP_005620.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P48029
Protein Weight	71kDa
Shipping	Ice bag