

Anti-Slc6A9 Polyclonal Antibody For Wb And Elisa - P48067

Item Number: CM0019464



Introduction

Rabbit polyclonal antibody targeting human SLC6A9 glycine transporter (GLYT1). Suitable for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Predicted molecular weight: 78 kDa. Ideal for neuroscience studies.

[Learn More](#)

Parameter	Value
Product Name	SLC6A9 Rabbit pAb
Remarks/Alias	GLYT1; GCENSG; SLC6A9
Species	Human
Gene ID (Human)	6536
Gene ID	6536
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human SLC6A9 (NP_964012.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48067
Protein Weight	78kDa
Shipping	Ice bag