

Anti-Slc7A11/Xct [Kd Validated] Rabbit Monoclonal Antibody For Wb - Q9Upy5

Item Number: CM0008179



Introduction

Validated rabbit monoclonal antibody targeting SLC7A11/xCT (cystine/glutamate transporter). Suitable for WB, IHC-P, ELISA with human, mouse, rat reactivity. Recombinant protein immunogen. 55 kDa. Ideal for ferroptosis and redox research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] SLC7A11/xCT Rabbit mAb
Remarks/Alias	xCT; CCBR1; SLC7A11/xCT
Species	Human
GeneID	23657
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UPY5
Protein Weight	55 kDa
Shipping	Ice bag