

Anti-Slc8A2 (Ncx2) Rabbit Polyclonal Antibody For Wb, Elisa - Q9Upr5

Item Number: CM0014580

Introduction

Rabbit polyclonal antibody against human SLC8A2 (NCX2), validated for WB and ELISA. Targets the sodium/calcium exchanger 2, a key regulator of neuronal Ca²⁺ homeostasis, synaptic plasticity, and urinary ion excretion.

[Learn More](#)

Parameter	Specification
Product Name	SLC8A2 Rabbit pAb
Remarks/Alias	NCX2; SLC8A2
Species	Human
GeneID (Human)	6543
GeneID	6543
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 360-450 of human SLC8A2 (NP_055878.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UPR5
Protein Weight	100kDa
Shipping	Ice bag