

Anti-Slc9A3R1 Rabbit Polyclonal Antibody For Wb, Elisa - O14745

Item Number: CM0013648



Introduction

Anti-SLC9A3R1 rabbit polyclonal antibody validated for WB and ELISA, targeting human, mouse, and rat NHERF1/EBP50 scaffold protein. Ideal for studying ion transport, Wnt signaling, and membrane-cytoskeleton interactions.

[Learn More](#)

Parameter	Value
Product Name	SLC9A3R1 Rabbit pAb
Remarks/Alias	EBP50; NHERF; NHE-RF; NHERF-1; NPHLOP2; SLC9A3R1
Species	Human
Gene ID (Human)	9368
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 189-358 of human SLC9A3R1 (NP_004243.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14745
Protein Weight	39kDa
Shipping	Ice bag