


## Anti-Slc9A3R2 Polyclonal Antibody For Wb, Elisa - Q15599

Item Number: CM0021505



### Introduction

Rabbit polyclonal antibody targeting SLC9A3R2 (NHERF2), validated for WB and ELISA. This scaffold protein links plasma membrane proteins to the actin cytoskeleton and regulates surface expression. Ideal for cell signaling research.

[Learn More](#)

| Parameter        | Value                                                                                                                   |
|------------------|-------------------------------------------------------------------------------------------------------------------------|
| Product Name     | SLC9A3R2 Rabbit pAb                                                                                                     |
| Remarks/Alias    | SIP1; OCT52; SIP-1; TKA-1; E3KARP; NHE3RF2; NHERF-2; SLC9A3R2                                                           |
| Species          | Human                                                                                                                   |
| Gene ID (Human)  | 9351                                                                                                                    |
| Gene ID          | 9351                                                                                                                    |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 180-300 of human SLC9A3R2 (NP_001123484). |
| Source           | Rabbit                                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                                   |
| Application      | WB, ELISA                                                                                                               |
| Cross Reactivity | Human                                                                                                                   |
| SWISS            | Q15599                                                                                                                  |
| Protein Weight   | 37kDa                                                                                                                   |
| Shipping         | Ice bag                                                                                                                 |