

Anti-Slc9A6 Rabbit Polyclonal Antibody For Wb And Elisa - Q92581

Item Number: CM0027447

Introduction

SLC9A6 rabbit polyclonal antibody for WB and ELISA, cross-reacts with mouse. Detects sodium/hydrogen exchanger 6, involved in endosomal pH regulation and neurodevelopment. Ideal for neuroscience and ion transport research.

[Learn More](#)

Parameter	Value
Product Name	SLC9A6 Rabbit pAb
Remarks/Alias	MRSA; NHE6; MRXSCH; SLC9A6
Species	Human
GeneID (Human)	10479
GeneID	10479
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 542-701 of human SLC9A6 (NP_001036002.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q92581
Protein Weight	78kDa
Shipping	Ice bag