

Anti-Slco1A2 Polyclonal Antibody For Wb, Elisa - P46721

Item Number: CM0028617

Introduction

Rabbit polyclonal antibody against human SLCO1A2 (OATP1A2), validated for WB and ELISA, cross-reactive with rat. For studying organic anion transport, drug uptake, and bile acid homeostasis.

[Learn More](#)

Parameter	Value
Product Name	SLCO1A2 Rabbit pAb
Remarks/Alias	OATP; OATP-A; OATP1A2; SLC21A3; SLCO1A2
Species	Human
Gene ID (Human)	6579
Gene ID	6579
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 400-500 of human SLCO1A2 (NP_066580.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P46721
Protein Weight	74kDa
Shipping	Ice bag