


Anti-Slco1B1 Polyclonal Antibody For Wb, Elisa - Q9Y6L6

Item Number: CM0024583



Introduction

SLCO1B1 Rabbit pAb, validated for WB and ELISA, detects human SLCO1B1 (OATP-C) with cross-reactivity to mouse and rat. Explore hepatic organic anion transport and statin pharmacology with this 76 kDa transporter antibody.

[Learn More](#)

Parameter	Description
Product Name	SLCO1B1 Rabbit pAb
Remarks/Alias	LST1; HBLRR; LST-1; OATP2; OATPC; OATP-C; OATP1B1; SLC21A6; SLCO1B1
Species	Human
GeneID (human)	10599
GeneID	10599
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 597-691 of human SLCO1B1 (NP_006437.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6L6
Protein Weight	76kDa
Shipping	Ice bag