

## Anti-Slco2B1 Rabbit Polyclonal Antibody For Wb/Elisa - 094956

Item Number: CM0014529

### Introduction

Rabbit polyclonal antibody against human SLCO2B1 (OATPB/OATP2B1), suitable for WB and ELISA. Detects ~77kDa organic anion transporter in human and mouse; ideal for steroid/drug transport and hepatic/placental studies.

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | SLCO2B1 Rabbit pAb                                                                                                   |
| Remarks/Alias    | OATPB; OATP-B; OATP2B1; SLC21A9; SLCO2B1                                                                             |
| Species          | Human                                                                                                                |
| GeneID (Human)   | 11309                                                                                                                |
| GeneID           | 11309                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 130-210 of human SLCO2B1 (NP_009187.1) |
| Source           | Rabbit                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                |
| Application      | WB, ELISA                                                                                                            |
| Cross Reactivity | Human, Mouse                                                                                                         |
| SWISS            | O94956                                                                                                               |
| Protein Weight   | 77kDa                                                                                                                |
| Shipping         | Ice bag                                                                                                              |