

Anti-Slco3A1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Uig8

Item Number: CM0016075

Introduction

CamelBio's SLCO3A1 rabbit polyclonal antibody (OATP-D, OATP3A1) is validated for WB and ELISA. Detects human, mouse, and rat SLCO3A1. Ideal for organic anion transporter research and tissue distribution studies.

[Learn More](#)

Parameter	Value
Product Name	SLCO3A1 Rabbit pAb
Remarks/Alias	OATPD; OATP-D; OATP3A1; OATPRP3; OATP-RP3; SLC21A11; SLCO3A1
Species	Human
GeneID (Human)	28232
GeneID	28232
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UIG8
Protein Weight	77kDa
Shipping	Ice bag