

Anti-Slco6A1 Polyclonal Antibody For Wb, Elisa - Q86Ug4

Item Number: CM0015746

Introduction

Rabbit polyclonal antibody against human SLCO6A1 (OATP6A1, CT48), validated for WB and ELISA. Detects human, mouse, and rat. Ideal for cancer/testis antigen and organic anion transporter research.

[Learn More](#)

Parameter	Value
Product Name	SLCO6A1 Rabbit pAb
Remarks/Alias	GST; CT48; OATPY; OATP-I; OATP6A1; SLCO6A1
Species	Human
GeneID (Human)	133482
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 480-580 of human SLCO6A1 (NP_775759.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86UG4
Protein Weight	79kDa
Shipping	Ice bag