

Anti-Slfn11 Rabbit Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - Q7Z7L1

Item Number: CM0006498



Introduction

High-quality SLFN11 rabbit monoclonal antibody validated for WB, IF/ICC, IP, and ELISA. Targets human SLFN11 (103 kDa) for DNA damage, antiviral, and replication stress research.

[Learn More](#)

Parameter	Value
Product Name	SLFN11 Rabbit mAb
Remarks/Alias	SLFN8/9
Species	Human
Gene ID (Human)	91607
Gene ID	91607
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 132-328 of human SLFN11 (NP_689483.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human
SWISS	Q7Z7L1
Protein Weight	103 kDa
Shipping	Ice bag